

ANTIMICROBIAL INSTANT-READ THERMOMETERS

MODEL SPECIFICATIONS



3621

8018

9842

9847

9848

9998

Model 3621

- Range 0/220°F
- Accuracy +/- 2°F
- NSF approved
- 5" stem with 1" dial face
- Exclusive Safe-T-Guard case inhibits growth of bacteria
- Exclusive Safe-T-Set system with built-in pocket wrench

Model 8018

- Range 0/220°F
- Accuracy +/- 2°F
- NSF approved
- 5" stem with 1-3/4" dial face
- Exclusive Safe-T-Guard case inhibits growth of bacteria
- Exclusive Safe-T-Set system with built-in pocket wrench

Model 9842

- Range -40/450°F, -40/230°C
- Accuracy +/- 1°F in HACCP ranges
- NSF approved
- Exclusive Safe-T-Guard case inhibits growth of bacteria
- Field recalibratable, waterproof case
- On/Off switch with auto shut-off after 10 minutes
- Rapid response time
- F&C selectable

Model 9847

- Range -40/450°F, -40/230°C
- Accuracy +/- 1°F in HACCP ranges
- NSF approved
- Exclusive Safe-T-Guard case inhibits growth of bacteria
- Waterproof case
- Pen style design
- Auto shut-off after 10 minutes
- Rapid response time with hold feature
- F&C selectable

Model 9848

- Range -40/450°F, -40/230°C
- Accuracy +/- 1°F in HACCP ranges
- NSF approved
- Exclusive Safe-T-Guard case inhibits growth of bacteria
- Waterproof and Dishwasher safe
- Max registry and field re-calibration features
- Rapid response time with hold feature
- On/Off switch with auto shut-off after 10 minutes
- F&C selectable
- Pen Style

Model 9997

- Safe-T-Guard Sanitizing System, bulk pack of 4 Sanitizing Tubes per box

Model 9998

- Saf-T-Guard Sanitizing system with model 3621 Antimicrobial Thermometer

EXCLUSIVE SAFE-T-GUARD CASE INHIBITS GROWTH OF BACTERIA AND FUNGUS

The fact is, even after a thermometer is thoroughly cleaned and disinfected, there is the potential of cross contamination from the pocket case.

Taylor's exclusive Safe-T-Guard sleeve incorporates a non-toxic chemical compound that inhibits the growth of harmful bacteria, which minimizes cross contamination. It is important to thoroughly clean and disinfect your thermometer after each use. Use Taylor Model 9997 Sanitizing System.

For Customer Service: 800-225-4834